

**Accessories****Active SmartProbes**

- AI 401 x10 low-noise differential amplifier
- AI 402 x50 ultra-low-noise differential amplifier
- AI 405 x10 differential amplifier

**SmartProbe Adapter**

- AI 417 passive adapter for ECG/EMG/EEG 2 mm leads

**SmartProbe Kits**

- AI 490 connector kit
- AI 491 cable kit

# Smart eight-channel programmable signal conditioner.

**CyberAmp 380 Signal Conditioner**

This computer-controlled signal conditioner is used in conjunction with any laboratory A/D system. It provides 4-pole Bessel low-pass filtering, AC/DC coupling, notch filtering, variable gain, and baseline correction. There are eight independent channels. Each channel can accept a SmartProbe to adapt to a wide variety of electrodes or transducers. An audio monitor, chart recorder outputs and calibration signal outputs are included.

For PC users, a cable to attach to the computer's 9-pin serial port is included with each CyberAmp

## MDS Analytical Technologies/Axon CNS Ordering Information

Part No.	Product Description	Price
MultiClamp	MultiClamp 700B computer-controlled current & patch clamp amp	\$ 14,060
Axoclamp	Axoclamp 900A computer-controlled current & voltage clamp	\$ 9,394
Axopatch	Axopatch 200B-2 capacitor feedback patch clamp amp	\$ 12,035
SoftPanel	SoftPanel (USB) optional control panel	\$ 983
Digidata	Digidata 1440A data acquisition system	\$ 5,261
pCLAMP	pCLAMP 10 Standard electrophysiology software (Windows)	\$ 4,759
pCLAMP Upgd	pCLAMP 10 Upgrade available for previous versions of pCLAMP	\$ 433+
CyberAmp	CyberAmp 380-2 eight-channel signal conditioner	\$ 7,873
Axoprotector	Axoprotector 800A single-cell electroporator	\$ 6,524
Mo-1-CV-7B	Headstage CV-7B patch clamp (standard) for MultiClamp 700B	\$ 2,812
Mo-HL-U	Electrode holder for all Universal (U)-type headstages	\$ 169
Mo-HS-9A-X10U	HS-9A headstage for Axoclamp 900A (choose x0.1, x1, x10 U)	\$ 476
Mo-Surcharge	Molecular Devices \$100 surcharge for orders under \$500	\$ 100
	Complete Axon CNS cellular neuroscience product line avail.	

U.S./Canada prices shown. International prices add 15%. Email or visit web store for latest prices.

