

"The EasyCode software you provided is simple to use. We use a large screen Macintosh for creating and storing protocols. We especially like the ability to create new protocols and download them directly to the ValveBank, it really couldn't be easier."

Dr. Susan Abrahamson
University of California,
Berkeley


EasyCode Dialog Box

- Enter start and stop times accurate to 0.01 sec.
- Use all ValveBank commands including multi-channel looping and interface triggering

Free your computer during experiments

Transfer and run programs stored in the ValveBank, which can be disconnected from the computer – leaving it available for data acquisition.

Expand the computing power of your ValveBank system.



EasyCode® Software

Program your ValveBank® with a Macintosh or PC using easy "click-and-drag" time bars. Save and load unlimited programs to disk, print out program listings, then download your sequences into the ValveBank's memory with the included serial cable in one simple step.

- **Program up to sixteen channels of valves and digital outputs**
Open multiple experiment windows. List sequences to screen or printer. Copy and paste. Zoom in and out of your experiment.

EasyCode Software Ordering Information

Part No.	Product Description	Price
01-02-USB	EasyCode® Mac - ValveBank programming software Macintosh OS 6-9 & Classic with USB-to-serial adapter	\$ 274
01-03	EasyCode® PC - ValveBank programming software - Windows 3.1-95	\$ 195
01-03-2000	EasyCode® PC - ValveBank programming software - Windows 98-XP	\$ 249
USB	USB to serial converter - PC or Macintosh	\$ 79

U.S./Canada prices shown. International prices add 15%. Email or visit web store for latest prices.

AutoMate Scientific®
READY FOR RESEARCH.™

800.998.MATE | www.autom8.com | 650 University Ave #5, Berkeley, CA 94710 USA
tel 510.845.6283 | fax 510.665.3975 | e-mail info@autom8.com