

**General Specifications**

**Net Load Capacity**

150 lbs, 67 kg

**Facilities Required**

80 psi air or nitrogen

**Finish**

Electropolished Stainless Steel  
Sides and Base. Textured Powder  
Coat Carbon Steel Top Plate

**Dimensions**

18(w) x 20(d) x 3"(h)  
(458 x 508 x 76mm)

**Weight**

50 lbs (22kg)

Lightweight system with exceptional low frequency passive vibration isolation.



**66 SERIES Active TableTop™ CSP®**

High-Performance Vibration Isolation System with Compact Sub-Hertz Pendulum Technology

**Benefits**

- Self-leveling design with height control valves
- Electropolished stainless steel chassis with powder finish steel top plate
- Incorporates patented CSP isolator design
- Lightweight, low-profile

**Features**

- Exceptional passive vibration isolation, comparable to our full-size industry standard 63-500 Series Tables
- Lightweight, Compact Design (less than 50 lbs) is easily portable

**TMC Air Table Ordering Information**

Part No.	Product Description	Price
TM-66-501	TableTop CSP Isolation Systems	\$ 2,290
	Complete TMC product line available. Call or email for part numbers.	

U.S./Canada prices shown. International prices add 15%. Email or visit web store for latest prices.

